

GENESY 3000 Rehab



For your wellness, body and beauty care. Ideal for demanding sportspeople.

194 Programs

16 Beauty- 11 Rehabilitation – 33 Tens- Pain –
23 Microcurrents - 9 Denervated - 1 Ionophoresis –
3 Incontinence – 6 Interferentials 92 Electrotherapy Menu
(Ems - Russian)

NEW! **"EASY PROGRAM" ELECTROTHERAPY MENU:** a section where it is possible to rapidly set a program by setting the desired parameters (time, frequency, body area) and to execute specific treatments using EMS and Russian (Kotz) currents.

MULTI-USER: it permits to store up to 25 different users and to create a library of personalized programs depending on individual aims.

AUTOMATIC MODE (AUTO STIM®): it permits to execute a program in automatic mode with no need of intervening manually in the intensity regulation.

SYNCROSTIM: It is possible to increase the intensity of the 4 channels simultaneously using only one button .

RUN TIME FUNCTION: this function permits to modify the time, the frequency and the pulse amplitude during the treatment

PROGRAMMABLE: 15 free memories, only for EMS and TENS.

2+2 MODE: it permits to execute two different programs simultaneously.

FAVORITES: it permits to store up to 15 programs for each user

LAST 10: it stores the last 10 executed programs

WORKING TIME: It indicates the total usage time of the device

MY TRAINER: it shows the electrode position on the display

STIM LOCK



TECHNICAL FEATURES

Channels: 4 independent CHANNELS
Frequency: 0,3 – 150 Hz
Power: 120 mA per channel
Pulse amplitude: 40 – 450 µs
Wave form: Square, biphasic, symmetrical and compensated wave
Power supply: mains or battery
Device in II A class
Backlit display
In compliance with the directive 93/42/CEE-2007/47/CEE

EQUIPMENT

1 Carrying case
1 GLOBUS GENESY 3000 Rehab Stimulator
4 Cables for electrode connection
2 cables for microcurrent
4 self-adhesive square electrodes
4 self-adhesive rectangular electrodes
1 Charger
1 Operating manual

Genesy 3000 Rehab

Program list

REHABILITATION

Swollen ankles
Atrophy recovery
Hemiplegy-upper limbs
Hemiplegy-lower limbs
Recovery after ACL surgery
Functional recovery
Ankle re-education
Leg re-education
Muscle Spasms
Shoulder subluxation prevention
Vastus medialis reinforcement
Motor point pen
Agonist/Antagonist
Muscle reinforcement

PAIN

High-frequency antalgic tens
Knee pain
Endorphinic tens
Post-surgical pain
Chronic pain
Shoulder pain
Carpal tunnel
Trapezius pain
Knee osteoarthritis
Rotator cuff tendinitis
Low frequency antalgic tens
Muscle pain
Nerve compression
Spinal compression pain
Spinal osteoarthritis
Spinal osteoporosis
Hip osteoarthritis
Ankle osteoarthritis
Muscle tendon injury pain
Post-surgical pain (bones fractures and soft tissue injuries)
Sciatica
Chronic lumbago
Cervical pain
Acute pain inguinal hernia
Bursitis-tendinitis
Muscle injuries
Bone fractures
Whiplash
Epicondylitis

IONOPHORESIS

Ionophoresis

INCONTINENCE

Mixed incontinence
Stress incontinence
URGE incontinence

INTERFERENTIALS

Knee osteoarthritis
Acute lumbago
Cervical pain
Chronic lumbago
Frozen shoulder
Knee post-surgical pain

MICROCURRENTS

Epicondylitis
Scapulohumeral periarthritis
Muscle restoration
Contusion
Edema
Skin ulcer
Sciatica
Lumbago
Brachial neuralgia
Acute pain
Articular pain
Stiff neck
Whiplash
Cervical spondylosis
Shoulder sprain
Carpal tunnel
Knee sprain
Knee osteoarthritis
Ankle sprain
Achilles tendon inflammation
Patella tendon inflammation
Rotator cuff inflammation
Tendon inflammation

DENERVATED

Triangular 1
Triangular 2
Triangular 3
Trapezoidal 1
Trapezoidal 2
Trapezoidal 3
Rectangular 1
Rectangular 2
Rectangular 3

BEAUTY

Drainage
Lipolysis
Toning massage
Skin tone improvement
Connective massage
Post-pregnancy drainage
Post-pregnancy lipolysis
Post-pregnancy firming

ELECTROTHERAPY

The Electrotherapy program list includes 48 EMS and 44 Russian(KOTZ) parameters combinations