

ACTIVA 700



For your wellness, body and beauty care. Ideal for sportspeople.

271 Programs

13 Face-Skin G-Pulse © – **60** Beauty – **58** Fitness – **53** Sport
21 Special Sports – **12** Tens – **3** Rehabilitation – **3** Incontinence
12 Microcurrents – **36** "3S" Serial Sequential Stimulation

ADVANCED TECHNOLOGIES IN THE NEW 2013 RANGE

Globus presents an innovative technology intended to revolutionize the electrostimulation application in beauty, sport and medical fields: **3S (Serial Sequential Stimulation) technology.**

"3S" SERIAL SEQUENTIAL STIMULATION: The sequential contraction of the different muscular groups produces a deep pressure wave in the musculature involved causing the interstitial fluid drainage and enhancing the return of the venous blood to the heart.

FACE-SKIN TREATMENTS: ACTIVA 700 offers the new innovative and regenerative **G-PULSE** currents, proposed for Aesthetic and Beauty applications and able to regenerate the collagen production in the tissues.

G-TRODE HEAD: ACTIVA 700 is equipped with a special head that permits, using the G-PULSE programs, to obtain a more effective skin stimulation.

NEW TRENDY CARRYING CASE: ACTIVA 700 is supplied with a practical modern carrying case in aluminum and eco leather that can be carried in an easy way and that makes you feel always fashionable.

FACE AND BREAST AREA: ACTIVA 700 is equipped with a complete set of programs for face and breast toning

SPECIAL SPORTS: a special program list to stimulate in an effective and specific manner the muscular areas used in the more common sport activities.

MULTI-USER: it permits to store up to 10 different users and to create a library of personalized programs depending on individual aims.

AUTOMATIC MODE (AUTO STIM®): it permits to execute a program in automatic mode with no need of intervening manually in the intensity regulation.

SYNCROSTIM: It is possible to increase the intensity of the 4 channels simultaneously using only one button.

RUN TIME FUCTION: this function permits to modify the time, the frequency and the pulse amplitude during the treatment

PROGRAMMABLE: 10 free memories, only for EMS and TENS.

2+2 MODE: it permits to execute two different programs simultaneously.

FAVORITES: it permits to store up to 15 programs for each user

LAST 10: it stores the last 10 executed programs

WORKING TIME: It indicates the total usage time of the device

MY TRAINER: it shows the electrode position on the display

STIM LOCK



TECHNICAL FEATURES

Channels: 4 independent CHANNELS

Frequency: 0,3 – 150 Hz

Power: 120 mA per channel

Pulse amplitude: 40 – 450 µs

Wave form: Square, biphasic, symmetrical, compensated

Power supply: rechargeable battery

Back-lit display

Device in II A class

In compliance with the directive

93/42/CEE-2007/47/CEE

EQUIPMENT

1 carrying case

1 ACTIVA 700 electrostimulator

1 G-Trode head

4 electrostimulation cables and 2 microcurrent cables

4 self-adhesive square electrodes and 4 self-adhesive rectangular electrodes

1 Gel bottle

1 Charger

1 Operating manual

ACTIVA 700

Program list

Fitness

Firming
Bio-Pulse firming
Sculpting
Bio-Pulse sculpting
Toning
Mass building
Body sculpting
Definition
Jogging
Anaerobic fitness
Aerobic fitness
Cramp prevention

Beauty

Drainage
Bio-Pulse drainage
Lipolysis
Bio-Pulse relaxing massage
Toning massage
Energizing massage
Connective massage
Swollen arms
Face capillaries
Definition
Slim figure
Skin tone improvement
Post-pregnancy drainage
Post-pregnancy lipolysis
Post-pregnancy firming
Breast firming
Breast sculpting

Sport

Capillarization
Warm-up
Pre-competition warm-up
Active recovery
Maximum strength
Endurance strength
Explosive strength
Reactivity
Post-competition recovery
Hypertrophy

Rehabilitation

Quadriceps atrophy (with knee prosthesis)
Recovery after ACL surgery
Shoulder subluxation prevention

Pain-Tens

Conventional antalgic tens
Menstrual pain
Modulated antalgic tens
Knee pain
Endorphinic tens
Cronich pain
Shoulder pain (s. h. syndrome)
Muscle pain
Chronic lumbago
Cervical pain
Bursitis-tendinitis
Osteoarthritis

G-Pulse

Breast microlifting
Cleavage microlifting
Face microlifting
Skin elasticity
Cellulite
BioSkin collagen
Wrinkles
Eye wrinkles
Expression wrinkles
Stretch marks
Nail strenghtening
Bio peeling
Hematoma

Special Sports

Running
Skiing
Golf

"3s" Serial Sequential Stimulation

The "3S" program list includes 36 parameter combinations

Area	Name	Hz	Delay time
Up.Limb	SerSeqStim 0,5 s	30	0,5
	SerSeqStim 1 s		1
Lo.Limb	SerSeqStim 2 s	50	2
	SerSeqStim 3 s		3
Trunk	SerSeqStim 4 s		4
	SerSeqStim serial		11
Total	36 Programs		

The "3S" programs are characterized by a delay activation of the channels 3 and 4 compared with the channels 1 and 2. The Serial Sequential Stimulation permits to stimulate the musculature in kinetic chain thanks to the differentiated activation times of the muscular groups involved.

"Delay time" refers to the delay seconds that the next pulse needs to start.

Microcurrents

Epicondylitis
Scapulohumeral periarthritis
Contusion
Sciatica
Lumbago
Acute pain
Articular pain
Stiff neck
Whiplash
Shoulder sprain
Knee sprain
Tendon inflammation

Incontinence

Mixed incontinence
Stress incontinence
URGE incontinence